






ORTHOTECH
Innovative Orthopedic Implant Manufacturer

Multi-Fix[®]
Proximal Femur Nail
System

Multi-Fix[®] System
Development
Cooperation
Technology

Anatomically Designed Orthopedic Trauma Implants
Innovative Technology Promising High Quality Products





Instruments and Implants

Plate Ordering Information

골 내로 진입 하기 용이한 형상, 환자 다양한 신체 크기 및 구조에 대응 가능하며 Target device 체결 상태에서 Lag screw를 고정할 수 있는 구조 Lag screw 의 슬라이딩이 가능한 구조이며 Lag screw 통하여 파절된 bones 접합시키기 위해 본 시멘트 투입 Proximal femur nail End cap 다양한 길이 구성과 Distal locking screw 체결 시, Cortical bone에 체결되는 부분의 고정력을 증대시켜 일정기간 견고한 고정 상태를 유지하도록 하여 골절환자의 자연 치유를 유도

■ Multi-Fix® PF-Nail / 120°

Order Code	Description	Color
N5001	PF-Nail 9 × 170 (120°)	Non Anodizing
N5002	PF-Nail 9 × 200 (120°)	
N5003	PF-Nail 9 × 230 (120°)	
N5004	PF-Nail 10 × 170 (120°)	Purple
N5005	PF-Nail 10 × 200 (120°)	
N5006	PF-Nail 10 × 230 (120°)	
N5007	PF-Nail 11 × 170 (120°)	Blue
N5008	PF-Nail 11 × 200 (120°)	
N5009	PF-Nail 11 × 230 (120°)	
N5010	PF-Nail 12 × 170 (120°)	Yellow
N5011	PF-Nail 12 × 200 (120°)	
N5012	PF-Nail 12 × 230 (120°)	



■ Multi-Fix® PF-Nail / 125°

Order Code	Description	Color
N5031	PF-Nail 9 × 170 (125°)	Non Anodizing
N5032	PF-Nail 9 × 200 (125°)	
N5033	PF-Nail 9 × 230 (125°)	
N5034	PF-Nail 10 × 170 (125°)	Purple
N5035	PF-Nail 10 × 200 (125°)	
N5036	PF-Nail 10 × 230 (125°)	
N5037	PF-Nail 11 × 170 (125°)	Blue
N5038	PF-Nail 11 × 200 (125°)	
N5039	PF-Nail 11 × 230 (125°)	
N5040	PF-Nail 12 × 170 (125°)	Yellow
N5041	PF-Nail 12 × 200 (125°)	
N5042	PF-Nail 12 × 230 (125°)	



■ Multi-Fix® PF-Nail / 130°

Order Code	Description	Color
N5061	PF-Nail 9 × 170 (130°)	Non Anodizing
N5062	PF-Nail 9 × 200 (130°)	
N5063	PF-Nail 9 × 230 (130°)	
N5064	PF-Nail 10 × 170 (130°)	Purple
N5065	PF-Nail 10 × 200 (130°)	
N5066	PF-Nail 10 × 230 (130°)	
N5067	PF-Nail 11 × 170 (130°)	Blue
N5068	PF-Nail 11 × 200 (130°)	
N5069	PF-Nail 11 × 230 (130°)	
N5070	PF-Nail 12 × 170 (130°)	Yellow
N5071	PF-Nail 12 × 200 (130°)	
N5072	PF-Nail 12 × 230 (130°)	



Instruments and Implants

Screw Ordering Information



■ Multi-Fix® PFN-Lag Screw

Order Code	Description	Color
N5401	PFN-Lag Screw 70	Yellow
N5402	PFN-Lag Screw 75	
N5403	PFN-Lag Screw 80	
N5404	PFN-Lag Screw 85	
N5405	PFN-Lag Screw 90	
N5406	PFN-Lag Screw 95	
N5407	PFN-Lag Screw 100	
N5408	PFN-Lag Screw 105	
N5409	PFN-Lag Screw 110	



■ Multi-Fix® PFN-Lag Blade

Order Code	Description	Color
N5421	PFN-Lag Blade 70	Blue
N5422	PFN-Lag Blade 75	
N5423	PFN-Lag Blade 80	
N5424	PFN-Lag Blade 85	
N5425	PFN-Lag Blade 90	
N5426	PFN-Lag Blade 95	
N5427	PFN-Lag Blade 100	
N5428	PFN-Lag Blade 105	
N5429	PFN-Lag Blade 110	



■ Multi-Fix® PFN-Lag screw/blade Fix module

Order Code	Description	Color
N5491	PFN-Lag screw/blade Fix module	Non Anodizing



■ Multi-Fix® PFN-Inner Screw

Order Code	Description	Color
N5481	PFN-Inner Screw 9 x 21	Yellow

Instruments and Implants

Screw Ordering Information

■ Multi-Fix® PFN-End Cap

Order Code	Description	Color
N5471	PFN-End Cap 12.5 × 11	Yellow
N5472	PFN-End Cap 16.0 × 16	
N5473	PFN-End Cap 16.0 × 21	
N5474	PFN-End Cap 16.0 × 26	






■ Multi-Fix® PFN-Distal Locking Screw

Order Code	Description	Color
N5441	PFN-Distal Locking Screw 5.0 × 20	Blue
N5442	PFN-Distal Locking Screw 5.0 × 22	
N5443	PFN-Distal Locking Screw 5.0 × 24	
N5444	PFN-Distal Locking Screw 5.0 × 26	
N5445	PFN-Distal Locking Screw 5.0 × 28	
N5446	PFN-Distal Locking Screw 5.0 × 30	
N5447	PFN-Distal Locking Screw 5.0 × 32	
N5448	PFN-Distal Locking Screw 5.0 × 34	
N5449	PFN-Distal Locking Screw 5.0 × 36	
N5450	PFN-Distal Locking Screw 5.0 × 38	
N5451	PFN-Distal Locking Screw 5.0 × 40	
N5452	PFN-Distal Locking Screw 5.0 × 42	
N5453	PFN-Distal Locking Screw 5.0 × 44	
N5454	PFN-Distal Locking Screw 5.0 × 46	
N5455	PFN-Distal Locking Screw 5.0 × 48	
N5456	PFN-Distal Locking Screw 5.0 × 50	
N5457	PFN-Distal Locking Screw 5.0 × 52	
N5458	PFN-Distal Locking Screw 5.0 × 54	
N5459	PFN-Distal Locking Screw 5.0 × 56	
N5460	PFN-Distal Locking Screw 5.0 × 58	
N5461	PFN-Distal Locking Screw 5.0 × 60	
N5462	PFN-Distal Locking Screw 5.0 × 64	
N5463	PFN-Distal Locking Screw 5.0 × 68	
N5464	PFN-Distal Locking Screw 5.0 × 72	
N5465	PFN-Distal Locking Screw 5.0 × 76	
N5466	PFN-Distal Locking Screw 5.0 × 80	
N5467	PFN-Distal Locking Screw 5.0 × 85	
N5468	PFN-Distal Locking Screw 5.0 × 90	
N5469	PFN-Distal Locking Screw 5.0 × 95	
N5470	PFN-Distal Locking Screw 5.0 × 100	



Instruments and Implants

Instruments Ordering Information



■ PFN System Instruments Set

Order Code	Description
IN51-01	Ruler
IN51-02	Guide Wire 3.0
IN51-03	Entry Drill Guide Sleeve
IN51-04	Entry Trocar Sleeve
IN51-05	T-Handle
IN51-06	Entry Reamer
IN51-07	Awl
IN51-08	Ball Tip Guide
IN51-09	Flexible Reamer 9
IN51-10	Flexible Reamer 10
IN51-11	Flexible Reamer 11
IN51-12	Flexible Reamer 12
IN51-13	Flexible Reamer 13
IN51-14	Connecting Screw
IN51-15	Connecting Screw Driver
IN51-16	Insertion Handle
IN51-17	Impacter
IN51-18	Spanner
IN51-19	Hammer
IN51-20	Targeting Arm 120°
IN51-21	Targeting Arm 125°
IN51-22	Targeting Arm 130°
IN51-23	Targeting Sleeve
IN51-24	Compression Nut
IN51-25	Guide Sleeve
IN51-26	Lag Screw-Trocar
IN51-27	Anti Rotation Guide Sleeve
IN51-28	Lag Screw Depth Gauge
IN51-29	Lag Screw Starter Drill
IN51-30	Lag Screw Reamer
IN51-31	Stopper
IN51-32	Lag Screw Driver
IN51-33	Flexible Driver
IN51-34	Distal Locking Screw Sleeve
IN51-35	Distal Drill Sleeve
IN51-36	Distal Locking Trocar
IN51-37	Drill Bit
IN51-38	Distal Screw Depth Gauge
IN51-39	Distal Locking Screw Driver
IN51-40	Connecting Screw Remove Driver
IN51-41	Lag Screw Extractor
IN51-42	Extractor Hammer
IN51-43	Remove Connecting Screw
IN51-44	Nail Extractor
IN51-45	L-Wrench Clean
IN51-46	Wire
IN51-47	Aiming Device For Lag Screw



78 Yuram-ro, Dong-gu, Daegu, 41059, Rep.of Korea
Tel +82-53-314-7016 Fax +82-53-314-7017

Web <http://www.orthotech.co.kr>
E_mail orthotech@orthotech.co.kr