



ORTHOTECH
Innovative Orthopedic Implant Manufacturer

**Multi-Fix[®]
Small System
Plate&Screw**

Multi-Fix[®] System
Development
Cooperation
Technology

Anatomically Designed Orthopedic Trauma Implants
Innovative Technology Promising High Quality Products


Distal Radius Locking Plate

MFDS. Liscence# 13-1611

Material Pure Titanium Gr.2

EDI Code C5475080 ::: Reconstruction Straight Type

Order Distal Radlus Narrow Plate
Distal Radlus Broad Plate



Description	Color	Spec	Cat.No
Distal Radius Small Plate	Blue	1Hole Left	P2421
		2Hole Left	P2423
		3Hole Left	P2425
		1Hole Right	P2422
		2Hole Right	P2424
		3Hole Right	P2426
Distal Radius Narrow Plate	Yellow	1Hole Left	P2411
		2Hole Left	P2413
		3Hole Left	P2415
		1Hole Right	P2412
		2Hole Right	P2414
		3Hole Right	P2416
Distal Radius Broad Plate	Green	1Hole Left	P2401
		2Hole Left	P2403
		3Hole Left	P2405
		1Hole Right	P2402
		2Hole Right	P2404
		3Hole Right	P2406

MFDS. Liscence# 13-1611
17-388

Material Titanium 6AL-4V-ELI

EDI Code C6411228 ::: Locking Cortical Screw - Can. Type
C6402180 ::: Locking Cortical Screw
C6402080 ::: Cortical Screw

Spec	Lock. Cannulated	Lock. Cortical	Cortical
Ø2.5 x 12	S3142	S3122	S3107
Ø2.5 x 14	S3143	S3123	S3109
Ø2.5 x 16	S3144	S3124	S3111
Ø2.5 x 18	S3145	S3125	S3112
Ø2.5 x 20	S3146	S3126	S3113
Ø2.5 x 22	S3147	S3127	S3114
Ø2.5 x 24	S3148	S3128	S3115
Ø2.5 x 26	S3149	S3129	S3116
Ø2.5 x 28	S3150	S3130	S3117
Ø2.5 x 30	S3151	S3131	S3118
Color	Blue	Yellow	Grey



Instruments

Cat.No	Description
IP20-24	Depth Gauge 60mm
IP20-16	Locking Driver 2.0 Hex.
IP20-20	Small Plate Driver Ø2.5mm - ST Pick Up
IP20-03	Drill Bit 2.0 (2.5 Cortical Screw)
IP20-11	Drill Guide - Long 2.0
IP20-14	Double Drill Sleeve 2.0x2.7
IT20-04	Small Plate System Case Radius

Proximal Humerus Locking Plate


MFDS. Liscence# 13-1611

EDI Code C5475080

Material Pure Titanium Gr.2


Proximal Humerus Locking Plate
Dynamic washer

Order




Proximal Humerus Broad Plate



Color	Spec	Cat.No
	2Hole	P2101
	3Hole	P2102
	4Hole	P2103
	5Hole	P2104

Proximal Humerus Narrow Plate



Color	Spec	Cat.No
	2Hole	P2111
	3Hole	P2112
	4Hole	P2113
	5Hole	P2114

Clavicle Locking Plate

MFDS. Liscence# 17-387


Material Pure Titanium Gr.2

EDI Code

C5478128 ::: Clavicle Shaft Plate

C5478228 ::: Clavicle Hook Plate

C5478328 ::: Clavicle End Plate



Clavicle Hook Plate



Clavicle End Plate



Clavicle Shaft Plate



Color	Hole	Height	Left	Right
Blue	6Hole	12mm	P2001	P2002
		15mm	P2003	P2004
		18mm	P2005	P2006
		12mm	P2007	P2008
	7Hole	15mm	P2009	P2010
		18mm	P2011	P2012
		12mm	P2013	P2014
		8Hole	15mm	P2015
	18mm	P2017	P2018	

Color	Hole	Left	Right
Grey	6Hole	P2021	P2022
	7Hole	P2023	P2024
	8Hole	P2025	P2026

Color	Hole	Left	Right
Yellow	6Hole A	P2041	P2042
	6Hole B	P2043	P2044
	7Hole A	P2045	P2046
	7Hole B	P2047	P2048
	8Hole	P2049	P2050
	9Hole	P2051	P2052
	10Hole	P2053	P2054


Distal Fibular Locking Plate

MFDS. Liscence# 13-1611

EDI Code C5430080

Material Pure Titanium Gr.2

Distal Fibular Plate - Broad Type



Color	Spec	Cat.No
Green	7Hole	P4601
	8Hole	P4602
	9Hole	P4603
	10Hole	P4604
	11Hole	P4605
	12Hole	P4606
	13Hole	P4607

Distal Fibular Plate - Narrow Type



Color	Spec	Cat.No
Yellow	7Hole	P4611
	8Hole	P4612
	9Hole	P4613
	10Hole	P4614
	11Hole	P4615
	12Hole	P4616
	13Hole	P4617

Distal Fibular Plate - Hook Type



Color	Spec	Cat.No
Blue	6Hole	P4621
	7Hole	P4622
	8Hole	P4623
	9Hole	P4624



Use With
 Ø3.5mm Cortical Screw
 Ø3.5mm Locking Screw
 Ø4.0mm Locking Cortical Screw - Cannulated Type
 Ø4.0mm Cancellous Screw



Calcaneal Locking Plate



MFDS. Liscence# 13-1611

Material Pure Titanium Gr.2

EDI Code C5467480

Color	Spec	Cat.No
Grey	Small (12H)	P4701
	Medium (16H)	P4702
	Large (19H)	P4703

Small Plate System

MFDS. Licence# 06-346

Material Pure Titanium Gr.2

EDI Code

C5463080 ::: Small Compression
C5461080 ::: One-Third Locking Plate
C5469080 ::: Reconstruction Locking Plate

Small Plate - Compression Type



Small Plate - One Third Tubular Type



Small Plate - Reconstruction Straight Type




Description	Color	Spec	Cat. No
Small Compression	[Grey Bar]	4 Hole	F301
		5 Hole	F302
		6 Hole	F303
		7 Hole	F304
		8 Hole	F305
One-Third Locking Plate	[Grey Bar]	4 Hole	F311
		5 Hole	F312
		6 Hole	F313
		7 Hole	F314
		8 Hole	F315
Reconstruction Locking Plate	[Grey Bar]	5 Hole	F322
		6 Hole	F323
		7 Hole	F324
		8 Hole	F325



Instruments

Cat.No	Description
IP20-04	Drill Bit 2.0 (2.5 Cortical Screw)
IP20-06	Drill Bit 3.5 (4.0 Cancellous Screw)
IP20-08	Taps 3.5
IP20-09	Taps 4.0
IP20-29	Locking Driver 2.5Hex(Tip)
IP20-15	Double Drill Sleeve 3.5x2.5
IP20-10	T-Handle Top
IP20-26	Small Plate Bender Left
IP20-27	Small Plate Bender Right
IT20-01	Small Locking Plate Case
IP20-24	Depth Gauge 60mm
IP20-16	Locking Driver 2.0 Hex
IP20-20	Small Plate Driver Ø2.5mm - ST Pick Up (IP20-16)
IP20-17	Locking Driver 2.5 Hex
IP20-21	Small Plate Driver- ST Pick Up (IP20-17)
C590	Small Locking Plate Screw Case
C590-1	Small Locking Plate Box
IP20-03	Drill Bit 2.0 (2.5 Cortical Screw)
IP20-05	Drill Bit 2.7 (3.5 Cortical Screw)
IP20-11	Drill Guide - Long 2.0
IP20-13	Drill Guide - Long 2.7
IP20-14	Double Drill Sleeve 2.0 x 2.7
IS00-10	Torque Driver 1.0Nm
IS00-15	Torque Driver 1.5Nm
IS01-20	Torque Driver Tip 2.0 Hex
IS01-25	Torque Driver Tip 2.5 Hex




Anatomical Locking Plate System Screw

MFDS. Liscence# 13-1611 ::: Locking Cortical Screw
 17-388 ::: Locking Cortical Screw-Cannulated Type
 06-345 ::: Cortical / Cancellous Screw

EDI Code C6402080 ::: Cortical Screw
 C6405080 ::: Cancellous Screw
 C6402180 ::: Locking Cortical Screw
 C6411228 ::: Locking Cortical Screw-Cannulated Type
 D1201080 ::: Washer





Material Titanium 6AL-4V-ELI ::: Screw



Length	Ø3.5 Cortical Self	Ø 3.5 Locking Cortical Screw	Ø 4.0 Locking Cortical Screw Cannulated Type	Ø 4.0 Cancellous Screw Full Type	Ø 4.0 Cancellous Screw Half
12mm	S3502	S3522	S3702	S3732	S3752
14mm	S3503	S3523	S3703	S3733	S3753
16mm	S3504	S3524	S3704	S3734	S3754
18mm	S3505	S3525	S3705	S3735	S3755
20mm	S3506	S3526	S3706	S3736	S3756
22mm	S3507	S3527	S3707	S3737	S3757
24mm	S3508	S3528	S3708	S3738	S3758
26mm	S3509	S3529	S3709	S3739	S3759
28mm	S3510	S3530	S3710	S3740	S3760
30mm	S3511	S3531	S3711	S3741	S3761
32mm	-	S3532	S3712	-	-
34mm	-	S3533	S3713	-	-
35mm	S3514	-	-	S3742	S3762
36mm	-	S3534	S3714	-	-
38mm	-	S3535	S3715	-	-
40mm	S3517	S3536	S3716	S3743	S3763
45mm	S3518	S3537	S3717	S3744	S3764
50mm	S3519	S3538	S3718	S3745	S3765
Color	Grey	Blue	Yellow	Grey	Grey

KIPO Multi-Fix® Dual Pitch System

Intellectual Property / No. 10-2014-0113695





78 Yuram-ro, Dong-gu, Daegu, 41059, Rep.of Korea
Tel +82-53-314-7016 Fax +82-53-314-7017

Web <http://www.orthotech.co.kr>
E_mail orthotech@orthotech.co.kr

ORTHO TECH Information In Compliance With The
Korea Food and Drug Regulations Place a Notice That The Product Produced.

OTBR-54 2022-03-04 Rev1