




ORTHOTECH
Innovative Orthopedic Implant Manufacturer

**Multi-Fix[®]
HTO
Plate System**

Multi-Fix[®] System
Development
Cooperation
Technology

Anatomically Designed Orthopedic Trauma Implants
Innovative Technology Promising High Quality Products

HTO T Plate & BARO-FIX™ HTO Plate



Characteristic Features

01. 해부학적인 디자인으로 Bone과의 높은 밀착성.

02. Angular Stability로 뛰어난 안전성.

03. Block 1홀 + 3홀 함께 사용 가능
사용자에게 더 많은 옵션 제공.

01. 오른쪽 & 왼쪽 맞춤 PLATE 제공.

02. Plate 양쪽 홈 기능으로 MCL, Hamstring
골부착부 골막 등 재봉합할시 높은 만족도.

Multi-Fix HTO Plate Instrument

Material

Pure Titanium (ASTM F67) ::: Plate / TI-6AL-4V-ELI ALLOY (ASTM F136) ::: Screw / TI PURE GRADE2 (ASTM F67) ::: Wedge Plate

HTOT Locking Plate



BARO-FIX™ Plate



Right

Left

HTO Plate Set



Locking Screw



Cortical Screw



Locking Screw
Can Type



Wedge Spiked Block Bolt



Wedge Spiked Block



1Hole

3Hole

Plate List

Color	Cat.No	Hole
Yellow	P5015	T 3Hole
	P5016	T 4Hole
Blue	P50L03	BARO LEFT 3HOLE
	P50L04	BARO LEFT 4HOLE
	P50R03	BARO RIGHT 3HOLE
	P50R04	BARO RIGHT 4HOLE

Wedge Spiked Block Chart

Color	Size	1Hole	3Hole	
Grey	5T (4.5 - 5.5)	P5026	P50B05	
	6T (5.5 - 6.5)	P5027	P50B06	
	7T (6.5 - 7.5)	P5028	P50B07	
	8T (7.5 - 8.5)	P5029	P50B08	
	9T (8.5 - 9.5)	P5030	P50B09	
	10T (9.5 - 10.5)	P5031	P50B10	
	11T (10.5 - 11.5)	P5032	P50B11	
	12T (11.5 - 12.5)	P5033	P50B12	
	13T (12.5 - 13.5)	P5034	P50B13	
	14T (13.5 - 14.5)	P5035	P50B14	
	15T (14.5 - 15.5)	P5036	P50B15	
	Block Connect Bolt	P5025		

Screw Chart

Locking Screw		Cortical Screw		Locking Screw Can Type	
Cat.No	Length	Cat.No	Length	Cat.No	Length
SL50020T	Ø5.0x20	S4503	Ø4.5x20 ●	S4803H	Ø5.0x20 ●
SL50022T	Ø5.0x22	S4504	Ø4.5x22 ●	S4804H	Ø5.0x22 ●
SL50024T	Ø5.0x24	S4505	Ø4.5x24 ●	S4805H	Ø5.0x24 ●
SL50026T	Ø5.0x26	S4506	Ø4.5x26	S4806H	Ø5.0x26
SL50028T	Ø5.0x28	S4507	Ø4.5x28	S4807H	Ø5.0x28
SL50030T	Ø5.0x30	S4508	Ø4.5x30	S4808H	Ø5.0x30
SL50032T	Ø5.0x32	S4509	Ø4.5x32	S4809H	Ø5.0x32
SL50034T	Ø5.0x34	S4510	Ø4.5x34	S4810H	Ø5.0x34
SL50036T	Ø5.0x36	S4511	Ø4.5x36	S4811H	Ø5.0x36
SL50038T	Ø5.0x38	S4512	Ø4.5x38	S4812H	Ø5.0x38
SL50040T	Ø5.0x40	S4513	Ø4.5x40	S4813H	Ø5.0x40
SL50042T	Ø5.0x42	S4514	Ø4.5x42	S4814H	Ø5.0x42
SL50044T	Ø5.0x44	S4515	Ø4.5x44	S4815H	Ø5.0x44
SL50046T	Ø5.0x46	S4516	Ø4.5x46	S4816H	Ø5.0x46
SL50048T	Ø5.0x48	S4517	Ø4.5x48	S4817H	Ø5.0x48
SL50050T	Ø5.0x50	S4518	Ø4.5x50	S4818H	Ø5.0x50
SL50055T	Ø5.0x55	S4519	Ø4.5x55	S4819H	Ø5.0x55
SL50060T	Ø5.0x60	S4520	Ø4.5x60	S4820H	Ø5.0x60
SL50065T	Ø5.0x65	S4521	Ø4.5x65	S4821H	Ø5.0x65
SL50070T	Ø5.0x70	S4522	Ø4.5x70	S4822H	Ø5.0x70

Instrument

Cat.No	Description
IS08-01	Screw Driver Handle
IG42-3346	Double Drill Sleeve 4.5 / 3.7
IG42-20	Wire Sleeve 2.0
IG42-18	Guide Pin 1.8
IG42-02	Can. Depth Gauge
C208	Hammer
IP50-04	Depth Gauge
C570	Z Retractor
IP50-05	Space Retractor
IP50-06	Osteotomy 22mm
IP50-07	Osteotomy 15mm
IP50-08	Osteotomy Impactor
IP50-09	Osteotomy Angular
IP50-10	Sliding Hammer
IP50-12	Insert Holder 4.0
IP50-14	HTO Space
IG42-43	Drill Sleeve 4.3
IG42-38	Drill Sleeve 3.8
IS05-T25	Screw Driver Shaft-T25
IS03-H35	Screw Driver Shaft-Hex3.5
IS04-H35	Can. Screw Driver Shaft - Hex3.5
ID42-42A	Drill Bit ASS'Y-Can. Type 4.2
IS03-H25	Screw Driver Shaft Hex 2.5
ID42-37A	Drill Bit 3.7
IS42-06	Pick Up

75 ~ 100 mm (5mm Unit) ●

● Order

Multi-Fix HTO Plate Instrument

C208 Hammer



ISO8-01 Screw Driver Handle



IP50-04 Z Retractor



IP50-05 Space Retractor



IP50-14 HTO Space



IG42-3346 Double Drill Sleeve 4.5 / 3.7



IP50-10 Sliding Hammer



IP50-12 Insert Holder 4.0



IG42-02 Can. Depth Gauge



IP50-06 Osteotomy 22mm



IP50-07 Osteotomy 15mm



IP50-08 Osteotomy Impacter



IP50-09 Osteotomy Angular



IS42-06 Pick Up



IG42-20 Wire Sleeve 2.0



IG42-43 Drill Sleeve 4.3



IG42-38 Drill Sleeve 3.8



IG42-18 Guide Pin 1.8



ID42-42A Drill Bit ASS' Y-Can, Type 4.2



ID42-37A Drill Bit 3.7



IS05-T25 Screw Driver Shaft T25



IS03-H35 Screw Driver Shaft Hex, 3.5



IS04-H35 Can, Screw Driver Shaft Hex, 3.5



IS03-H25 Screw Driver Shaft Hex, 2.5



OTBR-131987K 2022-02-22 Rev4