




ORTHOTECH

# <Multi-Fix><sup>®</sup> Nail System

*Tibial Type*  
*Femoral Type*




Nail System / Material : TI-6AL-4V-ELI ALLOY (ASTM F136)

### Introduction

Intramedullary fixation is very appropriate for the long bone diaphyseal fractures including tibia and femur. The implant is made from medical grade Titanium alloy supporting its safety and stability.

## <Multi-Fix> Nail System - Tibial Type



#### Diameter Ø8 (Unreaming)

- N7501 Tibial Type Ø8 x 270mm
- N7502 Tibial Type Ø8 x 285mm
- N7503 Tibial Type Ø8 x 300mm
- N7504 Tibial Type Ø8 x 315mm
- N7505 Tibial Type Ø8 x 330mm

#### Diameter Ø9 (Unreaming)

- N7511 Tibial Type Ø9 x 270mm
- N7512 Tibial Type Ø9 x 285mm
- N7513 Tibial Type Ø9 x 300mm
- N7514 Tibial Type Ø9 x 315mm
- N7515 Tibial Type Ø9 x 330mm

#### Diameter Ø10 (Reaming)

- N7521 Tibial Type Ø10 x 270mm
- N7522 Tibial Type Ø10 x 285mm
- N7523 Tibial Type Ø10 x 300mm
- N7524 Tibial Type Ø10 x 315mm
- N7525 Tibial Type Ø10 x 330mm

#### Diameter Ø11 (Reaming)

- N7531 Tibial Type Ø11 x 270mm
- N7532 Tibial Type Ø11 x 285mm
- N7533 Tibial Type Ø11 x 300mm
- N7534 Tibial Type Ø11 x 315mm
- N7535 Tibial Type Ø11 x 330mm

#### Diameter Ø12 (Reaming)

- N7514 Tibial Type Ø12 x 270mm
- N7542 Tibial Type Ø12 x 285mm
- N7543 Tibial Type Ø12 x 300mm
- N7544 Tibial Type Ø12 x 315mm
- N7545 Tibial Type Ø12 x 330mm

#### Diameter Ø13 (Reaming)

- N7551 Tibial Type Ø13 x 270mm
- N7552 Tibial Type Ø13 x 285mm
- N7553 Tibial Type Ø13 x 300mm
- N7554 Tibial Type Ø13 x 315mm
- N7555 Tibial Type Ø13 x 330mm



Tibial Nail X-ray View



## <Multi-Fix> Nail System - Femoral Type



Femoral Nail X-ray View

### Diameter Ø10

N5501	Femoral Type Ø10 x 320mm
N5502	Femoral Type Ø10 x 340mm
N5503	Femoral Type Ø10 x 360mm
N5504	Femoral Type Ø10 x 380mm
N5505	Femoral Type Ø10 x 400mm
N5506	Femoral Type Ø10 x 420mm
N5507	Femoral Type Ø10 x 440mm

### Diameter Ø11


N5511	Femoral Type Ø11 x 320mm
N5512	Femoral Type Ø11 x 340mm
N5513	Femoral Type Ø11 x 360mm
N5514	Femoral Type Ø11 x 380mm
N5515	Femoral Type Ø11 x 400mm
N5516	Femoral Type Ø11 x 420mm
N5517	Femoral Type Ø11 x 440mm

### Diameter Ø12

N5521	Femoral Type Ø12 x 320mm
N5522	Femoral Type Ø12 x 340mm
N5523	Femoral Type Ø12 x 360mm
N5524	Femoral Type Ø12 x 380mm
N5525	Femoral Type Ø12 x 400mm
N5526	Femoral Type Ø12 x 420mm
N5527	Femoral Type Ø12 x 440mm

### Diameter Ø13

N5531	Femoral Type Ø13 x 320mm
N5532	Femoral Type Ø13 x 340mm
N5533	Femoral Type Ø13 x 360mm
N5534	Femoral Type Ø13 x 380mm
N5535	Femoral Type Ø13 x 400mm
N5536	Femoral Type Ø13 x 420mm
N5537	Femoral Type Ø13 x 440mm



## <Multi-Fix> Nail System - Screw & End Cap



### Locking Screw Ø 5.0

N0411	Ø 5.0 x 25mm
N0412	Ø 5.0 x 30mm
N0413	Ø 5.0 x 35mm
N0414	Ø 5.0 x 40mm
N0415	Ø 5.0 x 45mm
N0416	Ø 5.0 x 50mm
N0417	Ø 5.0 x 55mm
N0418	Ø 5.0 x 60mm
N0419	Ø 5.0 x 65mm
N0420	Ø 5.0 x 70mm
N0421	Ø 5.0 x 75mm
N0422	Ø 5.0 x 80mm

### Locking Screw Ø 6.0


N0431	Ø 6.0 x 25mm
N0432	Ø 6.0 x 30mm
N0433	Ø 6.0 x 35mm
N0434	Ø 6.0 x 40mm
N0435	Ø 6.0 x 45mm
N0436	Ø 6.0 x 50mm
N0437	Ø 6.0 x 55mm
N0438	Ø 6.0 x 60mm
N0439	Ø 6.0 x 65mm
N0440	Ø 6.0 x 70mm



### End Cap

N0301	End Cap L - 0
N0302	End Cap L - 5
N0303	End Cap L - 10

## <Multi-Fix> Nail System - Instruments



- This is a medical device  
 - Recycling of single-use medical device is prohibited  
 - User-sterile product

IN75-02	Pocket
IN75-04	Awl - Curved
IN75-05	Skin Protector
IN75-06	Guide Pin Holder
IN75-07	Target Device for Tibia
IN75-08	Supine Driver
IN75-10	Extractor Large Bolt
IN75-11	Wrench Spanner
IN75-12	Slide Hammer
IN75-13	Nail impactor
IN75-14	Drill Guide 4.5
IN75-15	Drill Guide 5.0
IN75-16	Drill Sleeve
IN75-17	Sleeve Trocar
IN75-18	Drill Bit 4.5
IN75-19	Drill Bit 5.0
IN75-21	Depth Gauge
IN75-22	Awl
IN75-23	Driver
IN75-24	Universal Wrench
IN75-25	Ball Tip Guide
IN75-26	Straight Guide
IN75-28	Flexible Reamer Shaft
IN75-29	Reamer Tip 9.0
IN75-30	Reamer Tip 9.5
IN75-31	Reamer Tip 10.0
IN75-32	Reamer Tip 10.5
IN75-33	Reamer Tip 11.0
IN75-34	Reamer Tip 11.5
IN75-35	Reamer Tip 12.0
IN75-36	Reamer Tip 12.5
IN75-37	Reamer Tip 13.0
IN75-38	Reamer Tip 13.5
IN75-36	Nail Instrument Case
IN55-01	Pocket
IN55-02	Target Device
IN55-03	Drill Guide 6.0
IN55-04	Drill Bit 6.0
IN55-05	Reamer Tip 14.0
IN55-06	Reamer Tip 14.5

**ORTHOTECH**

- 78 Yuram-ro, Dong-gu, Daegu, 41059, Rep.of Korea
- Tel : +82-53-314-7016 • Fax : +82-53-314-7017
- Web : <http://www.orthotech.co.kr>
- E\_mail : [orthotech@orthotech.co.kr](mailto:orthotech@orthotech.co.kr)

OTBR-130608E-F 2022-03-03 rev3

