

**ORTHOTECH**  
Innovative Orthopedic Implant Manufacturer

**Multi-Fix<sup>®</sup>  
Art Band Pin**

Multi-Fix<sup>®</sup> System  
Development  
Cooperation  
Technology

Anatomically Designed Orthopedic Trauma Implants  
Innovative Technology Promising High Quality Products

# Art Band Pin

MFDS. Licence# 18-772

EDI Code BC4100AS :: Art Band Pin

Material STAINLESS STEEL - ASTM F138

Order


A pin used in Olecranon, Malleolar and Patella fractures, the pin is developed and made to be more convenient than the pre-existing tension pins.

## Indications

- Olecranon Osteotomies
- Transverse Olecranon Fractures
- Medial Malleolar Fractures
- Lateral Malleolar Fractures
- Patella Fractures



## Art Band Pin



Cat.No	A	B	C	D
E144	120	70	1.6	1.7
E148	140	70	2.0	2.1

## Snapper



Cat.No	A	B	C
IS80-02	112	8	2.2

## Case



OTBR-18772E 2022-03-03 Rev0